

NEXT GENERATION RT & DNA POLYMERASES

High Fidelity EasyScript™/plus RT

0.5kb 1.0kb 1.5kb 2.6kb 3.3kb 4.8kb 6.5kb 8.7kb



EasyScript™ Elongation Ability


PCR amplification using Human cDNA synthesized with EasyScript™, followed by electrophoresis on 1% agarose gel.

- The Best Quality
- The Most Reliable Products
- The Highest Fidelity
- The Best Prices

Cat. No.	Size	Price
G231	25 rxns	\$35.00
G177	25 rxns	\$55.00

High Fidelity Bestaq™ DNA Polymerases

1.5kb 2.6kb 3.3kb 4.8kb 6.5kb 8.7kb 12.4kb 15.6kb



Bestaq™ Specifically Amplifies Genomic Template DNA up to 15.6 kb

PCR amplification with Bestaq™ of various targets, ranging from 1.5 kb to 15.6 kb, from genomic DNA, followed by electrophoresis on a 1% agarose gel.

- The Best Quality
- Elongate at 1.0 kb/sec.
- The Highest Fidelity
(> 50X Taq Polymerase)

Cat. No.	Size	Price
G456	100U (50µl)	\$70.00
G457	400U (200µl)	\$235.00



1-866-757-2414 (Toll Free)
604-247-2416 (Local)



Order: order@abmGood.com
Website: www.abmGood.com



Applied Biological Materials Inc.
8 - 13520 Crestwood Place,
Richmond, BC, CANADA V6V 2G2