

QuickRun™ 15 Min. Running Buffer

Innovative Protein Electrophoresis Workflow Essentials






Fast forward your research & seeing sharper bands



Trick: by changing just one simple ingredient – the running buffer

- **Fast Speed:** 85% time saving from 90min. to 15 min. (@250V) compared to traditional buffer
- **High Performance:** Sharp and razor crisp protein bands as result of proprietary reduced heat generating buffer formulation that prevents protein modification or degradation
- **Compatibility:** Most precast gels* in market and homemade gels with all gel formulations of pH > 8,0

*exclude NuPAGE® Protein Gel. NuPAGE® is registered trade make of Life Technologies.



Order Information

- ACT-ID1591 QuickRun™ running buffer 10X stock solution, denaturing, with SDS
- ACT-ID1581 QuickRun™ running buffer 10X stock solution, native, without SDS



Email: sales@bio-rev.com
www: www.bio-rev.com
Phone:
 Singapore 65-62733022
 Malaysia 6012-5250822

ISO 9001:2008