


DirectLoad™ Precast Protein Gel

Innovative Protein Electrophoresis Workflow Essentials

Imagine directly load dirty crude sample onto gel and see sharp bands?

DirectLoad™ precast protein gel is a revolutionary high-performance polyacrylamide gel uniting unparalleled features and benefits in one product as result of Innovations of Protein Electrophoresis Proprietary Technology.

- **DirectLoad™: Only precast gel in market** for hassle-free crude sample (cell lysate and etc.) direct load while delivering razor sharp crisp bands. Eliminates sample purification step all together.
- **Simple Universal Gel for Gradient like separation: One fits all for protein 10-250kDa**, eliminates complication of different % gel or gradient gel.
- **Storage & Shelf Life: One year room temperature** or 18 months 4 degree shelf life.



Order Information

ACT-IS111	DirectLoad™ precast polyacrylamide for 10-250 kDa proteins, 10X10 cm 10 well
ACT-IS112	DirectLoad™ precast polyacrylamide for 10-250 kDa proteins, 10X10 cm 15 well
ACT-IS111s	DirectLoad™ precast polyacrylamide for 10-250 kDa proteins, 10X8.2 cm 10 well
ACT-IS112s	DirectLoad™ precast polyacrylamide for 10-250 kDa proteins, 10X8.2 cm 15 well



Email: sales@bio-rev.com
www.bio-rev.com
Phone:
Singapore 65-62733022
Malaysia 6012-5250822

ISO 9001:2008