



Saliva Field Kit

Catalog# 0706 RSID-Field Kit - 5 Test Packs/Kit

The new **R**apid **S**tain **I**dentification Saliva Field Kit (RSID®-Saliva Field Kit) is designed for fast, easy, and reliable detection of human saliva from evidence encountered at a crime scene. All of the materials needed to perform the stain analysis as well as to document and transport the results are included in the kit. RSID-Saliva Field Kit results correlate with the ability to obtain DNA profile data.

All Materials Needed for Stain Identification are Included

- Extraction/running buffer, water for swab wetting, swabs, scissors, transfer pipettes, tube rack, protocols, and RSID-Saliva strip tests
- Single use components eliminate cross contamination

Efficient Sample Analysis

- Stain sampling and extraction: 50 minutes.
- Strip test results: 10 minutes, Total time = 1hr

Record, Document, and Store

- Convenient form for reporting and photographing test results and evidence details
- Sealable envelope for evidence storage, transport, and chain of custody maintenance

Detect as Little as 1 μ l of Human Saliva

- RSID results correlate with the ability to obtain DNA profiles
- No cross reaction observed with human blood, semen, or urine
- No cross reaction observed with animal saliva

Immuno-Chromatographic Lateral Flow Strip Test

- Standard assay format
- Clear and easy to read results
- Extended shelf life

Complete Set of Protocols Included

Available Now!

- Stain sampling/cutting protocol
 - Forensic sample extraction method
-

4600 Roosevelt Rd., Ste. 201, Hillside, IL 60162
TEL. (866) 434-2400 - FAX (708) 978-5115
www.ifi-test.com/rsid.html